



Scanning Aperture Ion Beam Range Compensator and Intensity Modulator for External Radiation Therapy

[View U.S. Patent No. 7,977,648 in PDF format.](#)

WARF: P07279US

Inventors: Thomas (Rock) Mackie, David Westerly, Ryan Flynn

The Wisconsin Alumni Research Foundation (WARF) is seeking commercial partners interested in developing a technique for external beam ion radiation therapy that modulates an area beam by using a shutter system to increase treatment speed while maintaining high precision.

Overview

External beam radiation therapy is a tumor treatment technique that directs one or more high-energy radiation beams to the tumor. Advanced external beam radiation systems treat tumors with multiple X-ray fan beams that can be rotated around the patient. Each beam consists of individual “beamlets” whose intensity can be controlled to allow for the treatment of complex tumor shapes. However, x-ray energy still is deposited to all tissue along the entire path of the beam up to the exit point.

An alternate beam radiation therapy technique involves the use of ions such as protons to treat tumors. This technique dramatically reduces the radiation dose to healthy tissue when compared with X-ray radiation therapy using photons. The dose intensity from protons is not uniform along the beam path and rises to a “Bragg peak” near a point where the proton beam stops completely. Controlling the placement of the Bragg peak so it is located on the tumor reduces the amount of radiation delivered to the patient’s healthy tissue. Unlike X-ray radiation therapy, ions allow for separate control of the total dose of radiation (intensity) and distance the Bragg peak occurs (range).

Existing proton therapy techniques utilize one of two methods. The spread out Bragg peak (SOBP) approach spreads the ion beam out into an “area beam” using scattering foils and brass collimators to cover the entire tumor in one exposure. This method is fast, but is less precise and requires a special compensator to adjust for the tumor shape. The magnetic spot scanning (MSS) approach uses the original ion “pencil beam” and steers it around the tumor with successive exposures. Due to the multiple exposures the process is slow. The potential for high accuracy is negated due to patient movement, which risks missing tumor areas.

The Invention

UW–Madison researchers have developed a technique for external beam ion radiation therapy that modulates an area beam by the use of a shutter system to increase treatment speed while maintaining high precision. The shutters can control beam intensity independently to eliminate the need for a special compensator.

This technique uses an area beam that is occluded with a set of longitudinally opposed and latitudinally adjacent shutter pairs. Each shutter of each pair is controllable and can extend to different longitudinal distances. The shutters can control the range of the beam to

We use cookies on this site to enhance your experience and improve our marketing efforts. By continuing to browse without changing your browser settings to block or delete cookies, you agree to the storing of cookies and related technologies on your device. [See our privacy policy.](#)

Applications

OK



WARF
Wisconsin Alumni Research Foundation

| info@warf.org | 608.960.9850

- Ion beam radiation therapy for the treatment of tumors

Key Benefits

- Allows for independent control of intensity and range
- Increases accuracy of beam placement
- Increases treatment speed

Additional Information

Related Technologies

- [For a modulation system to be used in conjunction with this technology, see WARF reference number P07283US.](#)

Related Intellectual Property

- [View Divisional Patent in PDF format.](#)

Tech Fields

- [Radiation Therapy : External beam therapy.](#)

For current licensing status, please contact Jeanine Burmania at jeanine@warf.org or 608-960-9846

We use cookies on this site to enhance your experience and improve our marketing efforts. By continuing to browse without changing your browser settings to block or delete cookies, you agree to the storing of cookies and related technologies on your device. [See our privacy policy.](#)

OK



WARF
Wisconsin Alumni Research Foundation

| info@warf.org | 608.960.9850