



## Correction of Inverse Consistency and Transitivity Errors in Deformable Image Registration

[View U.S. Patent No. 8,526,692 in PDF format.](#)

**WARF: P110232US01**

Inventors: Wolfgang Tome, Edward Bender, Nicholas Hardcastle

**The Wisconsin Alumni Research Foundation (WARF) is seeking commercial partners interested in developing a corrective strategy to minimize discrepancies resulting from inverse-inconsistent and non-transitive deformation vector fields in radiotherapy dose accumulation.**

### Overview

In radiation therapy, particularly with 4-D planning and adaptive radiotherapy, multiple images almost always are acquired and utilized as part of the therapy. A treatment plan is defined on an image of the patient and used to align and control radiation beams during treatment. The dose received by each region of tissue during each treatment fraction is monitored. The registration of the treatment plan with the patient and the monitoring of radiation dose over time are complicated by changes in the patient and target over time. It is common to obtain new images of the patient at the time of each treatment fraction.

Deformable image registration (DIR) is necessary for accurate dose accumulation between multiple radiotherapy image sets. DIR algorithms can exhibit inverse and transitivity inconsistencies, which cause differences in accumulated dose from different image pathways. Differences in accumulated dose may impact clinical decisions. Deformation algorithms to reduce inverse consistency and transitivity errors exist; however, there is a need for post-processing techniques that minimize inconsistencies in accumulated doses.

### The Invention

UW–Madison researchers have developed a correction algorithm that forces inverse consistency and transitivity between all deformation vector fields used in dose accumulation. The algorithm provides flexible image deformation that is largely indifferent to the sequence of deforming step and does not require transformation through a particular image. These results are accomplished by ensuring individual deformation maps are inverse consistent, for example, by use of a correlation algorithm or an averaging process. Transitive deformation maps are produced by combining inverse-consistent deformation maps for combinations of different deformation pathways between the same endpoint images. The resulting system allows arbitrary deformation of image-type data among images with consistent results when the deformation endpoints are the same.

### Applications

- Dose-guided radiotherapy

### Key Benefits

- Allows flexibility in changing the registration of images from one to different images without hitting your browser cookies; you agree to the storing of cookies and related technologies on your device. [See our privacy policy.](#)
- Reduces the difference in accumulated dose between different image pathways
- Improves accuracy of delivered dose accumulation

OK



### Publications

- Bender E., Hardcastle N. and Tome W. 2012. On the Dosimetric Effect and Correction of Inverse Consistency and Transitivity Errors in Dose Accumulation. Med. Phys. 39, 272.

### Tech Fields

- [Radiation Therapy : Treatment planning](#)

For current licensing status, please contact Jeanine Burmania at [jeanine@warf.org](mailto:jeanine@warf.org) or 608-960-9846

We use cookies on this site to enhance your experience and improve our marketing efforts. By continuing to browse without changing your browser settings to block or delete cookies, you agree to the storing of cookies and related technologies on your device. [See our privacy policy.](#)

OK



**WARF**  
Wisconsin Alumni Research Foundation

| [info@warf.org](mailto:info@warf.org) | 608.960.9850