



Multiplexed Mass Spectrometry Quantification with Neutron-Encoded Mass Tags

[View U.S. Patent No. 9,366,678 in PDF format.](#)

WARF: P130032US01

Inventors: Joshua Coon, Alexander Hebert

The Wisconsin Alumni Research Foundation (WARF) is seeking commercial partners interested in developing methods that combine the best aspects of SILAC and isobaric tagging for more accurate proteomic quantification in a large number of samples.

Overview

In the past decade, cataloguing the thousands of proteins present in a cell using mass spectrometry (MS) has matured into a procedure both deeply vital and straightforward. Yet discovering how such molecules change in abundance under various circumstances remains a challenge.

A key step in most quantification strategies is to flag molecules of interest with stable isotopes that then can be detected and tallied by MS analysis. Two main approaches seek to accomplish this. SILAC (stable isotope labeling with amino acids in cell culture) incorporates 'light' or 'heavy' amino acids into proteins during cell growth. However, multiplexing (testing four or more samples simultaneously) is arduous with SILAC and redundant mass/charge peaks can lead to spectral complexity. Also, the technique is incompatible with mammalian biofluids and tissues. Isobaric tagging is a chemical labeling approach developed to overcome these issues, but interference from other molecules tends to degrade accuracy.

To empower modern proteomics and medicine, a new approach must support expanded, multiplexed analyses with heightened accuracy and reproducibility.

The Invention

UW-Madison researchers have developed customized tagging reagents for MS proteome quantification. Called neutron encoded mass tagging, or 'NeuCode,' the method exploits small mass differences, on the order of millidaltons (mDa), between heavy isotopes, such as nitrogen or carbon, which can be coupled with amino acids used during cell culture similarly to SILAC. Through instrument resolution settings, NeuCode labeling allows the user to control when multiplexed information is masked, limit sample complexity, and when it is revealed, produce quantitative data. Furthermore, these mDa mass differences can be incorporated into novel chemical tags that couple to peptides and enable analysis of biological samples from any source. Unlike isobaric tags, however, specialized reporter groups, linkers and charge sites are not required.

Applications

- Relative or absolute quantification of analytes, including amino acids, metabolites, peptides and proteins

Key Benefits

- Empowers high throughput analysis
- Useful for large number of samples at once
- Compatible with biological tissues and fluids

We use cookies on this site to enhance your experience and improve our marketing efforts. By continuing to browse without changing your browser settings to block or delete cookies, you agree to the storing of cookies and related technologies on your device. [See our privacy policy.](#)

OK



WARF
Wisconsin Alumni Research Foundation

| info@warf.org | 608.960.9850

- Alleviates spectral complexity problem
- No precursor interference caused by co-isolating impurities
- Does not require MS/MS scan
- Dissociation can be accomplished by any fragmentation method.

Stage of Development

An initial version of the tag has been tested and shown to enable 12- and 18-plex analysis at current MS resolving powers with the potential of 27-plex testing if resolution is doubled.

Additional Information

For More Information About the Inventors

- [Joshua Coon](#)

Related Technologies

- [WARF reference number P06069US describes ionizable isotopic labeling reagents for relative quantification by mass spectrometry.](#)

Tech Fields

- [Analytical Instrumentation, Methods & Materials : Mass spectrometry.](#)
- [Research Tools : Detection](#)

For current licensing status, please contact Jennifer Gottwald at jennifer@warf.org or 608-960-9854

We use cookies on this site to enhance your experience and improve our marketing efforts. By continuing to browse without changing your browser settings to block or delete cookies, you agree to the storing of cookies and related technologies on your device. [See our privacy policy.](#)

OK



WARF
Wisconsin Alumni Research Foundation

| info@warf.org | 608.960.9850