



Methods of Manufacturing and Using Beta-Peptide Lyotropic Liquid Crystals

[View U.S. Patent No. 8,021,570 in PDF format.](#)

WARF: P06005US

Inventors: Samuel Gellman, Nicholas Abbott, William Pomerantz

The Wisconsin Alumni Research Foundation (WARF) is seeking commercial partners interested in developing liquid crystals based on beta-peptide scaffolds.

Overview

Liquid crystals have been widely explored for applications such as display and sensing technologies. Lyotropic liquid crystals, which become more ordered as they become more concentrated within a solvent, are particularly useful for detecting the presence of biological targets, including proteins, viruses and microbes. Some alpha-peptides have been used to form lyotropic liquid crystals, but these peptides must be quite long, which limits their use.

The Invention

UW-Madison researchers have developed liquid crystals based on relatively short beta-peptide scaffolds that can be tailored to have unique properties for use in biological assays. Beta-peptides differ from conventional peptides in that certain beta-peptides form stable helices at short, oligomeric lengths, giving rise to robust, asymmetric structures. These helical beta-peptides can self-assemble to form lyotropic liquid crystals in aqueous environments.

To give the liquid crystals more desirable properties, a variety of functional groups can be added to the beta-peptides. The liquid crystals can then be used, for example, to detect an analyte, such as a protein, in a biological sample.

Applications

- Detection of proteins, cells or viruses

Key Benefits

- Modular synthesis of beta-peptides allows control over the location of functionalized side chains, including ionic groups that confer aqueous solubility and aromatic groups that influence liquid crystal properties.
- Oligomers of beta-amino acids may incorporate chemical groups and/or small molecules to further functionalize the resulting liquid crystal.
- Beta-peptides can be tailored to have specific viscosity and birefringence, important qualities for liquid crystals.
- Beta-peptides are resistant to proteolysis and generally non-toxic to cells, making them useful for biological applications, such as biomolecular sensing of proteins, cells or viruses in biological samples.

Additional Information

We use cookies on this site to enhance your experience and improve our marketing efforts. By continuing to browse without changing your browser settings to block or delete cookies, you agree to the storing of cookies and related technologies on your device. [See our privacy policy.](#)

For More Information About the Invention

- [Samuel Gellman](#)

OK



Tech Fields

- [Analytical Instrumentation, Methods & Materials : Sensors](#)
- [Research Tools : Detection](#)

For current licensing status, please contact Rafael Diaz at rdiaz@warf.org or 608-960-9847

We use cookies on this site to enhance your experience and improve our marketing efforts. By continuing to browse without changing your browser settings to block or delete cookies, you agree to the storing of cookies and related technologies on your device. [See our privacy policy.](#)

OK



WARF
Wisconsin Alumni Research Foundation

| info@warf.org | 608.960.9850