



Accurate Determination of Metabolite Concentrations in Cell Extracts or Biological Fluids

[View U.S. Patent No. 7,847,547 in PDF format.](#)

WARF: P07362US

Inventors: John Markley, Ian Lewis, William Westler, Qiu Cui, Gabriel Cornilescu, Marco Tonelli

The Wisconsin Alumni Research Foundation (WARF) is seeking commercial partners interested in developing a simple experimental protocol and device for identifying and measuring concentrations of metabolites in complex solutions.

Overview

Metabolomics, or the systematic study of all metabolites in an organism, is one of the fastest growing areas of scientific research. Nuclear magnetic resonance (NMR) spectroscopy is one tool used to analyze metabolites. Unlike other detection methods, the metabolites do not need to be separated for NMR, making the sample recoverable for further testing.

However, traditional one dimensional (1-D) NMR cannot identify compounds in complex mixtures of biological fluids and tissue extracts that may contain hundreds of overlapping peaks. This method is only accurate to within 16 percent in a simple mixture of known compounds.

Two dimensional (2-D) ^1H - ^{13}C NMR is a well-established technique for reducing peak overlap. But it is generally used for qualitative, rather than quantitative, analyses of metabolites. It is difficult to translate multidimensional peak intensities into metabolite concentrations, and 2-D spectra take more time to collect than 1-D spectra, making it impractical for metabolomics studies with hundreds of samples.

The Invention

UW-Madison researchers have developed a simple experimental protocol and device for identifying and measuring concentrations of metabolites in complex solutions from 2-D ^1H - ^{13}C NMR spectra. This protocol may lead to a faster, more consistent and accurate method for generating more complete metabolite concentration information.

The intensity peaks from the NMR data are used to identify metabolites in advance. Then a processor calculates the concentration of each metabolite compound, based on the concentration equation associated with the compound and the intensity value from the NMR data.

Applications

- Identifying and quantifying most NMR observable metabolites

Key Benefits

- Reduces acquisition times to as little as eight minutes per spectrum

We use cookies on this site to enhance your experience and improve our marketing efforts. By continuing to browse without changing your browser settings to block or delete cookies, you agree to the storing of cookies and related technologies on your device. [See our privacy policy.](#)

- Effective for identifying compounds and measuring molar concentrations of metabolites to within 2% percent, in contrast to traditional 1-D NMR

- Very close to a commercial product—requires little further development work



WARF
Wisconsin Alumni Research Foundation

| info@warf.org | 608.960.9850

Stage of Development

This protocol was successfully tested using a set of 30 synthetic mixtures containing varying concentrations of 26 small molecule metabolites that mimic extracts from *Arabidopsis thaliana*.

Additional Information

For More Information About the Inventors

- [John Markley](#)

Tech Fields

- [Research Tools : Genomics & proteomics](#)

For current licensing status, please contact Jennifer Gottwald at jennifer@warf.org or 608-960-9854

We use cookies on this site to enhance your experience and improve our marketing efforts. By continuing to browse without changing your browser settings to block or delete cookies, you agree to the storing of cookies and related technologies on your device. [See our privacy policy.](#)

OK



WARF
Wisconsin Alumni Research Foundation

| info@warf.org | 608.960.9850