



Microwave Cavity Detector for Highly Sensitive Mass Spectrometry Detection of Large Molecules

[View U.S. Patent No. 8,669,521 in PDF format.](#)

WARF: P100106US01

Inventors: Robert Blick, Jonghoo Park

The Wisconsin Alumni Research Foundation (WARF) is seeking commercial partners interested in developing an add-on device for mass spectrometers to allow highly sensitive time and mass measurements.

Overview

Mass spectrometers are analytic instruments that precisely measure the mass of molecules, and therefore may be used to determine composition and molecule identity based on mass data. Generally, the molecules are given an electrical charge (i.e., ionized) and then accelerated by an electrical field. The velocity of their acceleration will be proportional to the mass to charge ratio, and the molecular mass may be determined using a known charge and a velocity measurement. Velocity measurements are obtained using a variety of methods including use of a "sector" type analyzer or "time of flight" analyzer.

Mass spectrometry is applied increasingly to extremely large molecules such as proteins. The large mass of protein molecules decreases the sensitivity of a time of flight spectrometer because the velocity of the proteins is lower, and therefore the difference between velocities of similar masses is less. As a result, the time of flight measurement must be obtained with higher precision. Conventional spectrometer detectors are limited by measurement delays that can hide small mass differences between large molecules, which restricts their use in proteomics research that requires highly sensitive protein mass detection. Improvement in the mass resolving power of spectrometers is needed.

The Invention

UW-Madison researchers have developed a detector for mass spectrometry that employs a tuned microwave cavity to increase the sensitivity of time measurements. The detector can be used in mass spectrometers that provide a source of ionized molecules, which are accelerated in an electric field before reaching the detector. The design includes a cavity of conductive material providing an electromagnetically tuned cavity, which includes an opening positioned to receive the accelerated molecules. An antenna communicates with the cavity to receive an electrical signal caused by electromagnetic resonance of the cavity. Detection electronics receive this electrical signal, which marks the arrival of the ionized molecule in the cavity and allows for time of flight measurement. With proper configuration of the frequency of resonance in the cavity, its modes and its quality factor, time resolutions on the order of one nanosecond are possible.

Applications

- Mass spectrometry machines used in research such as proteomics that requires highly sensitive mass detection
- Detection and analysis of molecules with large molecular mass, such as proteins

We use cookies on this site to enhance your experience and improve our marketing efforts. By continuing to browse without changing your browser settings to block or delete cookies, you agree to the storing of cookies and related technologies on your device. [See our privacy policy.](#)

Key Benefits

OK



WARF
Wisconsin Alumni Research Foundation

| info@warf.org | 608.960.9850

- Improves sensitivity
- Improves temporal resolution
- Allows nanosecond time resolution

Additional Information

Related Technologies

- [For information about another tool to detect large molecular masses using mass spectrometry, see WARF reference number P05110US.](#)

Tech Fields

- [Analytical Instrumentation, Methods & Materials : Mass spectrometry.](#)
- [Research Tools : Detection](#)

For current licensing status, please contact Jennifer Gottwald at jennifer@warf.org or 608-960-9854

We use cookies on this site to enhance your experience and improve our marketing efforts. By continuing to browse without changing your browser settings to block or delete cookies, you agree to the storing of cookies and related technologies on your device. [See our privacy policy.](#)

OK



WARF
Wisconsin Alumni Research Foundation

| info@warf.org | 608.960.9850