



Wirelessly Tracking Interventional Medical Device with MRI System

[View U.S. Patent No. 8,805,475 in PDF format.](#)

WARF: P120271US01

Inventors: Krishna Kurpad, Madhav Venkateswaran

The Wisconsin Alumni Research Foundation (WARF) is seeking commercial partners interested in developing a method to transmit tracking and imaging signals from inside a patient using an internal coil powered by radiofrequency pulses.

Overview

Interventional medical procedures can be performed less invasively, and more safely, when guided by magnetic resonance imaging (MRI). Placing devices like catheters, guidewires and stents using MRI guidance has great clinical potential.

One particular challenge in the field, however, has been the development of safe and reliable means for tracking such devices as they are manipulated within vessels or organs. Every previous approach has drawbacks. Attaching and tracking X-ray markers requires radiation, for example, while incorporating paramagnetic material into devices degrades image quality.

Another approach is to connect the devices with special coils that give feedback in real time. But blurry imaging is a problem. Also, cables have to extend from the tip of the device, out of the patient, and into an external connection. A new MRI tracking system needs to be more workable, clear and accurate.

The Invention

UW–Madison researchers have developed a method for driving and amplifying MR signals using an imaging coil that is coupled to a medical device and inserted into a subject. The position-tracking coil uses the MR system's external radiofrequency pulses as a power source.

The coil switches between two duties: receiving and storing energy, and acquiring imaging data. Its signals can be wirelessly transmitted to the external MRI system using an amplifier and booster stages.

Applications

- MR-guided interventional procedures

Key Benefits

- No cumbersome power source
- Reduces the number of cables running to MRI scanner
- Less radiation exposure to patient than in X-ray applications
- Fewer adverse effects to patient from contrast agents than in X-ray applications

We use cookies on this site to enhance your experience and improve our marketing efforts. By continuing to browse without changing your browser settings to block or delete cookies, you agree to the storing of cookies and related technologies on your device. [See our privacy policy.](#)

Additional Information

OK



WARF
Wisconsin Alumni Research Foundation

| info@warf.org | 608.960.9850

Related Technologies

- [For more information about MR-guided devices for tracking and imaging in 3-D, see WARF reference number P06222US.](#)

Tech Fields

- [Medical Imaging : MRI](#)

For current licensing status, please contact Jeanine Burmania at jeanine@warf.org or 608-960-9846

We use cookies on this site to enhance your experience and improve our marketing efforts. By continuing to browse without changing your browser settings to block or delete cookies, you agree to the storing of cookies and related technologies on your device. [See our privacy policy.](#)

OK



WARF
Wisconsin Alumni Research Foundation

| info@warf.org | 608.960.9850