



Algorithms to Classify T Cell Activation by Autofluorescence Imaging

[View U.S. Patent No. 11,410,440 in PDF format.](#)

WARF: P190306US02

Inventors: Melissa Skala, Anthony Gitter, Zijie Wang, Alexandra Walsh

The Wisconsin Alumni Research Foundation (WARF) is seeking commercial partners interested in developing a method to non-invasively detect T cell activation by imaging NAD(P)H. These algorithms can be applied to NAD(P)H images taken with commercial imaging flow cytometers/sorters, and fluorescence microscopes.

Overview

The importance of T cells in immunotherapy is driving the development of technologies to better characterize T cells and improve therapeutic efficacy. One specific objective is assessing antigen-induced T cell activation because only functionally active T cells are capable of killing the desired targets. Autofluorescence imaging can assess functional activity of individual T cells in a non-destructive manner by detecting endogenous changes in metabolic co-enzymes such as NAD(P)H. However, recognizing robust patterns of T cell activity is computationally challenging in the absence of exogenous labels or information-rich autofluorescence lifetime measurements.

Non-invasive, contrast-agent free determination of T cell behavior is imperative for advancing studies of T cells *in vivo* and for quality control of clinical adoptive immune transfer procedures (e.g., CAR T cell therapies).

The Invention

Building on award-winning work, UW–Madison researchers have discovered that autofluorescence intensity images of NAD(P)H can accurately classify T cells as activated or not activated ('naïve' or 'quiescent'), and have developed algorithms to classify T cell activation based on the images. Specifically, adapting pre-trained convolutional neural networks (CNNs) for the T cell activity classification task, T cells can be classified with 92 percent accuracy. These pre-trained CNNs perform better than classification based on summary statistics (e.g., cell size) or CNNs trained on the autofluorescence images alone.

This invention provides a way to non-invasively detect T cell activation by imaging NAD(P)H intensity. These algorithms can be applied to NAD(P)H images taken with commercial imaging flow cytometers / sorters, and fluorescence microscopes. If increased accuracy of T cell activation is needed for a specific application, additional measurements of the other NAD(P)H and FAD fluorescence endpoints can be obtained and used for classification.

Applications

- Screening and sorting activated T cells

Key Benefits

- Algorithms can be applied to existing commercial systems

We use cookies on this site to enhance your experience and improve our marketing efforts. By continuing to browse without changing your browser settings to block or delete cookies, you agree to the storing of cookies and related technologies on your device. [See our privacy policy.](#)

- Method is non-invasive and contrast-agent free, hence can reduce reagent cost and variability
- Current methods to determine T cell activation require contrast agents and may require tissue/cell fixation

OK



WARF
Wisconsin Alumni Research Foundation

| info@warf.org | 608.960.9850

Stage of Development

The researchers have demonstrated that advanced machine learning models can accurately classify T cell activity from NAD(P)H intensity images and that those image-based signatures transfer across human donors.

Additional Information

For More Information About the Inventors

- [Melissa Skala](#)
- [Anthony Gitter](#)

Related Technologies

- [WARF reference number P180292US02 describes a method and device to detect T cell activation by detection of autofluorescence from free NAD\(P\)H, or NAD\(P\)H \$\alpha\$ 1. The device is also designed to allow sorting of T cells based on NAD\(P\)H \$\alpha\$ 1.](#)

Publications

- [Wang Z. J., Walsh A. J., Skala M. C. and Gitter A. 2020. Classifying T Cell Activity in Autofluorescence Intensity Images with Convolutional Neural Networks. Journal of biophotonics, 13\(3\), e201960050.](#)

Tech Fields

- [Analytical Instrumentation, Methods & Materials : Microscopy](#)
- [Medical Devices : Diagnostics & monitoring tools](#)

For current licensing status, please contact Jeanine Burmania at jeanine@warf.org or 608-960-9846

We use cookies on this site to enhance your experience and improve our marketing efforts. By continuing to browse without changing your browser settings to block or delete cookies, you agree to the storing of cookies and related technologies on your device. [See our privacy policy.](#)

OK



WARF
Wisconsin Alumni Research Foundation

| info@warf.org | 608.960.9850