



Imaging Spectrometer for Early Detection of Skin Cancer

[View U.S. Patent No. 8,315,692 in PDF format.](#)

WARF: P06477US

Inventors: Andrew Sheinis

The Wisconsin Alumni Research Foundation (WARF) is seeking commercial partners interested in developing a portable cancer detection system with improved imaging characteristics.

Overview

Skin cancer detection is typically identified via standard visual inspection. Skin patches with asymmetrical outlines, irregular borders, mottled color or a diameter greater than six millimeters are all early indicators of possible cancerous areas.

However, some skin cancers cannot be detected by visual inspection because they appear in light frequency bands other than the three bands visible to the human eye. For this reason, spectral information can be collected and analyzed to improve the early detection of skin cancer.

Currently, two techniques are used for spectrally imaging tissue. One method uses different frequencies of light, while the other uses a moveable scanning slit. In both cases, the imaging process is delayed, either by the need to switch between colors of light or by the need to scan a slit over the area of interest. This delay can cause the image reference points and light spectrum to become skewed, a problem that can be further exacerbated if movement occurs during the acquisition process. Additionally, the moving parts of current scanning devices make them cumbersome and unreliable, particularly for portable applications. New techniques for spectrally imaging tissue are needed.

The Invention

A UW-Madison researcher has developed a portable imaging spectrometer for the early detection of skin cancer. A handheld scanner uses light emitting diodes to illuminate a region of skin and the reflected light is collected by an objective lens. A micro-lens array then divides the region into smaller images that are processed to reveal their spectral content.

Because spectral and image data are acquired in one step, this new device provides two effective indicators to detect skin cancer. Physicians can evaluate the image data while the spectral data is compared to spectra of known cancerous or healthy regions.

Applications

- Skin cancer detection

Key Benefits

- Acquires high resolution spectral and image data in one step

We use cookies on this site to enhance your experience and improve our marketing efforts. By continuing to browse without changing your browser settings to block or delete cookies, you agree to the storing of cookies and related technologies on your device. [See our privacy policy.](#)

- Acquires a single image with multiple color channels
- Avoids image distortion

- Does not require moving components

OK



- Device is rugged and handheld.
- Samples all tissue within mapped region, so that small areas of cancer are not missed
- Accentuates regions with cancerous features
- Image is easily reviewed and evaluated by a physician.
- Image processing may take place within device or in a remote computer connected by a cable or wireless system.

Publications

- [Click here for a news release describing this technology.](#)

Tech Fields

- [Analytical Instrumentation, Methods & Materials : Spectroscopy](#)
- [Medical Devices : Diagnostics & monitoring tools](#)
- [Medical Imaging : Other diagnostic imaging](#)

For current licensing status, please contact Jeanine Burmania at jeanine@warf.org or 608-960-9846

We use cookies on this site to enhance your experience and improve our marketing efforts. By continuing to browse without changing your browser settings to block or delete cookies, you agree to the storing of cookies and related technologies on your device. [See our privacy policy.](#)

OK



WARF
Wisconsin Alumni Research Foundation

| info@warf.org | 608.960.9850